




# IMPACTO® FOOT PROTECTION

## Viscolas Modular Orthotic System

- > Prefabricated components are biomechanically engineered to provide effective treatment for a variety of common disorders
- > Components can be used in conjunction with Viscolas RELIEVER insoles or placed directly in the shoe
- > Viscolas® is the only visco-elastic polymer formulated for orthopedic application
- > Made of non-irritating Viscolas, resistant to fungal and bacterial growth and can be easily cleaned with soap and water
- > Natural tackiness of the components can be increased with an application of alcohol or decreased with talcum powder
- > Can be comfortably worn with dress shoes, sports shoes, work shoes and boots
- > Soft viscoelastic polymer won't "bottom out" with special surface tackiness adhesive to keep in place
- > Made in USA, PVC free



**METATARSAL PAD**  
Useful in conditions like metatarsalgia, painful plantar lesions, plantar fat pad atrophy, bursitis, capsulitis, periostitis, neurotrophic feet, Interdigital neuroma. Helps to reestablish normal weight bearing pattern.



**LONGITUDINAL ARCH PAD**  
Indications include plantar fasciitis, shin splints, heel pain, iliotibial band syndrome, foot and leg fatigue, medial ankle pain. Helps control excessive pronation or foot flattening.



**CORRECTIVE HEEL WEDGE**  
Indications include lateral instability. Bidirectional use  
Varus medial: reduces pronation. Valgus lateral decreases lateral instability.



**HEEL LIFT**  
Used for leg length discrepancy, Achilles tendonitis, equinus deformity.

PART No.	Item	Size	Dimensions (widest)
VI4942S	Metatarsal pad	Small	1 1/2" long x 1" wide
VI4943L	Metatarsal pad	Large	2 1/4" long x 1 5/8" wide
VI4944S	Arch pad	Small	2 1/2" long x 3/4" wide
VI4945M	Arch pad	Medium	3 1/4" long, 1 1/4" wide
VI4946L	Arch pad	Large	4" long, 1 1/2" wide
VI4999S	Heel wedge	Small 1/8"	3" long, 1 1/4" wide
VI4975L	Heel wedge	Large 1/8"	4" long, 1 3/4" wide
VI4976S	Heel wedge	Small 1/4"	3" long, 1 1/4" wide
VI4977L	Heel wedge	Large 1/4"	4" long, 1 3/4" wide
VI4985S	Heel lift	Small 1/8"	
VI4986L	Heel lift	Large 1/8"	
VI4997S	Heel lift	Small 1/4"	2 3/4" long, 2" wide
VI4998L	Heel lift	Large 1/4"	3 1/4" long, 2 1/4" wide

### PHYSICIANS COMPONENTS - MODULAR ORTHOTIC SYSTEM

Insoles and components are self-adhesive. Wide range of components allows flexibility to formulate immediate custom orthotic for patients requiring moderate adjustments. Combinations for therapeutic and preventive treatment.

TRUSTED PROTECTION

ISO 9001:2015



[IMPACTO.CA](http://IMPACTO.CA)

N. America 888 232 0031

Free phone UK 0800 0280 243

T: 613 966 0062 F: 613 966 0067

PO BOX 524, BELLEVILLE, ON K8N 5B2

[impacto@impacto.ca](mailto:impacto@impacto.ca)

Distributed by.

# IMPACTO

PROTECTIVE PRODUCTS