



IMPACTO® KNEE PROTECTION

Thermo Wrap Knee Support

Model TS209

- › Pull on style with knee patella cut-out
- › Soft lining wicks away moisture and helps prevent skin irritation
- › Reinforced stitching
- › Protects from Repetitive strain Injuries
- › Compression and support
- › Long term comfort
- › Lining captures body heat helping to increase circulation, allows skin ventilation
- › Durable multi-layered support



Part	Size	Circumference
TS20920	Small	12"-13 ³ / ₄ " 31-35cm
TS20930	Medium	14"-15 ³ / ₄ " 36-40cm
TS20940	Large	16"-17 ³ / ₄ " 41-45cm
TS20950	X-Large	18"-19 ³ / ₄ " 46-50cm
TS20960	XX-Large	20"-21 ³ / ₄ " 51-55cm

Sizing: Measure slightly bent underneath knee cap.
 Thermo Wraps can be used prior, during and after activities to help PREVENT RSI (Repetitive Strain Injuries). Fit should be firm but not tight. Wear lining directly against the skin. Do not use liniments, oils or ointments. Should an irritation or rash develop, discontinue use immediately and seek medical advice.

Washing: Machine wash in lukewarm water and mild detergent. Air dry. Do not bleach, do not dry by machine, do not dry clean. Wash daily if worn during heavy activity.



Sold per each (double for pair) 12 per Case

TRUSTED PROTECTION

ISO 9001:2015



IMPACTO.CA

N. AMERICA 888 232 0031
 FREE PHONE UK 0800 0280 243
 T. 613 966 0062 F. 613 966-0067
 PO BOX 524, BELLEVILLE, ON K8N 5B2
IMPACTO@IMPACTO.CA

Distributed by.

IMPACTO

PROTECTIVE PRODUCTS