

Manufacturer's Instructions and User Information

Name and Address of the manufacturer: Impacto Protective Products Inc.
PO Box 524, Belleville, ON K8N 5B2 Canada

Style Identification: VI4769, VI4770, VI4771, VI4772 Blue Dot

Description: **IMPACTO® Blue Dot Heel Cups**

Molded heel cups made of medical grade silicone designed with a softer "Blue Dot" in the medial side of the heel. This "Blue Dot" area is made of a softer silicone to help cushion the heel for greater comfort underfoot. The cupped design forces the fatty pad of the heel (the body's own natural defense against skeletal shock) downward for extra protection.



These garments bear CE marking to demonstrate compliance with EU 2016/425 Personal Protective Equipment Annex II Health & Safety Requirements.

EU Declaration of Conformity: <https://www.impacto.ca/products/datasheet-ce-conformity-statements/>

Fitting Instructions:

Remove any loose fitting inner insoles or debris from heel area of shoe. Place heel cup into the heel area of the shoe. Make sure the blue dot is on the inside of the heel.

Part#	USA Men	USA women	UK	Europe
VI4769S	2 – 4	4 – 6	1 – 3	33 - 36
VI4770M	5 – 7	7 - 9	4 – 6	37 – 40
VI4771L	8 – 10	10 - 12	7 – 9	41 – 44
VI4772XL	11 - 13		10 - 12	45 - 48

Weight (Lbs.): 0.30, Sold per Pair. Case: 50.

Washing Instructions:

Hand wash with mild soap and water. Air dry. Do not wash or dry by machine.

Limitations of Use:

It is recommended that, once inserted, these cushions remain in the shoes. Continual removal from one pair of shoes to another may cause deterioration.

After Care:

Follow the manufacturer's recommendations for washing and cleaning. To ensure optimum protection, regular comparisons are recommended against a new product. This product is warranted against defects in materials and workmanship for 12 months from purchase. There are no other warranties expressed or implied. The company will not accept liability for product, whose use and care advice has been ignored.

EU DECLARATION OF CONFORMITY

We,

Impacto Protective Products Inc.
PO Box 524
Belleville, ON K8N 5B2
Canada

Declare that the PPE described here after,

Blue dot, Heel Spur, Heel Cups, Reliever Insoles

These garments bear CE marking to demonstrate compliance with EU
2016/425 Personal Protective Equipment Annex II Health & Safety Requirements.

this declaration of conformity is self-certified and issued under the sole responsibility of the
manufacturer;

The technical Construction File is maintained at the following addresses:

1. Impacto Protective Products Inc.
PO Box 524
Belleville, ON K8N 5B2 Canada
Tel 0613 9660 062
Fax 0613 9660 067
UK free phone 0800 0280 243
Email impacto@impacto.ca
Contact Mr. Eric Lehtinen

Authorized Signatory



Date Oct 24, 24

Name:

Eric Lehtinen

Position:

President
Impacto Protective Products Inc.